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Glutaraldehyde (5,1%); HEMA (36,1%); Sodium Fluoride and Deionized Water.

#### **INDICATION:**

**G.H.F.** Desensitizer is indicated for treatment and prevention of hypersensitive dentin, mainly in the dental neck area; to be applied under restorations; post-surgery sensitivity; before impression (in patients with reported sensitivity to cold); after scaling and root planning; to reduce bacterial plaque; and to reduce cervical grinding progression caused by abrasion and erosion.

#### **TECHNICAL INFORMATION:**

**G.H.F.** Desensitizer is specially indicated for the treatment of hyperesthetic dentin because it recover and obliterate the dentin tubules, forming a barrier that inhibits the painful stimulus to pulp. This protection is intense on the first 6 months after application and a new application is required after this time. It has nano metric mineral loads which penetrate and reinforce dentin protecting it against cervical erosion. The result is instantaneous after the application and there is no restriction for ingestion of liquid and solid food. It's a non light cured material and forms no coating on dental surface. It's compatible with all adhesive system or restoration. It can be used as primer in cement or restorative system that needs this agent. The fluoride presents in its composition is an additional agent in dentin desensitizing.

#### **USE INSTRUCTIONS:**

##### **I - For sensitivity:**

1. Do the prophylaxis with pumice powder and water;
2. Isolate the region where **G.H.F.** Desensitizer will be placed and also protect the gingival tissue with solid Vaseline;
3. Dry using absorbing paper or cotton;
4. Place one drop of **G.H.F.** Desensitizer in a glass mixing pad and apply on all surface using a disposable brush;
5. Wait 30 seconds and dry smoothly with an air jet
6. Rinse with water.

##### **II - Applications under Restorations:**

1. Complete insulation (rubber dam) to protect mucosa;
2. Etching according to manufacturer's recommendation;
3. Wash all area with water and remove dentin humidity with a light airflow;
4. Pour one drop of **G.H.F.** Desensitizer in a glass mixing pad and apply on all surface using a disposable brush;
5. Wait during 30 seconds and dry it with a light airflow;
6. Apply Bonding (a.e. Master Bond) according to manufacturer's instructions and finish restoration.

#### **SPECIAL CARE:**

**G.H.F.** Desensitizer should not be in contact with the skin and mucosa; In case of eyes contact, wash