USE INSTRUCTIONS:

1 - Definitive esthetic restorations:

1. Do the crown polishing with pumice powder and water;
2. Select the shade with the tooth still moist;
3. Do the cavitory prepare in accordance with the chosen technique
4. Proceed the absolute isolation;
5. Do the etching (a.e. Attaque Gel), for 30 seconds at the enamel and for 15 seconds at the dentin;
6. Rinse thoroughly with water during 30 seconds;
7. Dry but do not dehydrate the dentin. It is important to dry well the enamel but do not dehydrate the dentin.
8. Apply the bonding agent (a.e. Master Bond) as indicated by the manufacturer's instructions;
9. Place **MASTER FILL** in thin layers (1 to 2mm) light curing each layer for 40 ± 10 seconds ;
10. Do the finishing ;
11. Check occlusion;
12. Do the polishing using abrasive discs and silicon abrasive tips;
13. It is recommended the use of a surface sealant for restoration (a.e. Bioforty) according to manufacturer's instructions.

II - Prosthesis Repair and Indirect Restorations:

1. To improve the adhesion, use the porcelain conditioner according to the manufacturer's instructions;
2. Apply the Silane on the broken porcelain and wait for 1 minute, do a new application;
3. Apply bonding system (a.e. Master Bond) according to the manufacturer's instruction;
4. Use **MASTER FILL** filling the cavity in thin layers (thickness of min 1 max 2 mm);
5. Do the finishing and polishing.

Observation: If necessary, disguise the metal.

SPECIAL CARES:

- The dentin-pulp system must be protected with an appropriate base before the etching;
- Eugenol based liners affect negatively in the light curing;
- The use of cationic collutorium, plaque indicators and clorexidine may interfere at the coloration of the material.
- To obtain a satisfactory polymerization, adequate equipments are required and the distance for the polymerization of the surface must be inferior of 5 mm. Follow rigorously the light exposure time, as well as the thickness of the layer; Non polymerized composites may cause irritation to the dentin-pulp system, changes in the shades and also reduce the mechanical properties.
- Absolute isolation is recommended.
- The application of surface sealant is recommended (a.e. Bioforty), in order to improve the life time of the restoration.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

STORAGE:

It is recommended to keep material inside refrigerator but product must be at ambient temperature before the use.

After using keep the syringes tightly closed to avoid exposure to the light.

Do not store near products containing Eugenol. Eugenol interferes in the polymerization.

CONTRAINDICATION:

MASTER FILL should not be used in patients is know to have sensibility to any of the ingredients listed.

WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.

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MASTER FILL

Light Curing composite for Restoration

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1x3g syringe Ataque Gel (37% Orthophosphoric Acid);

1x5mL flask Master Bond DE (Light Curing Single Component Adhesive for Enamel and Dentin);

20 Disposable Brushes;

1 Application Tip.

Reposition:

1x4g syringe **MASTER FILL**

Available in shades: A1-A2-A3-A3,5-A4-B0,5-B1-B2-B3-C1-C2-C3-C4-D3-INCISAL-OA2-OA3-OB2

Reg. ANVISA: 10298550040

COMPOSITION:

Bisphenol A glicidilmethacrylate; Urethane Ethildimethacrylate; Inorganic Filler; Pigments and Catalyst.

INDICATION:

It is a micro hybrid resin indicated for anterior and posterior definitive dental restorations in enamel and dentin; splinting of mobile teeth; prosthesis repair; indirect restorations.

TECHNICAL INFORMATION:

MASTER FILL is a high density composite resin for use in both anterior and posterior tooth restorations on enamel and dentin.

Its formulation contains inorganic particles varying between 0,04 and 2,2 µm.

It presents a total inorganic filler volume of 79%.

MASTER FILL presents fluorescence properties, excellent finishing, polishing and mechanics properties, allowing its use at anterior restorations, veneers and occlusal surfaces. Setting time of 40 ± 10 seconds, if using a 400 MW curing light (The indicated light curing range is 16 Joules).